

ENHANCE

Enhance Your System
with these Special Offers!



Connect Premium

NOW INCLUDED with our most popular system plans, for the ultimate protection and discounts for yourself and your equipment! Take advantage of this exclusive new program that provides you with an extra layer of protection AND can **save you up to \$300!**

- **Protect yourself** from costly replacements of lost or damaged equipment through our Connect Premium Program
- **Up to 50% off** select features, accessories and services
- **25% off** Fall Detect Pendant service
- **One FREE** standard Help Button

Only \$6/month additional
(if not included in your plan)



Fall Button™

Add the Fall Button for \$11/month in addition to your plan. Our unique Fall Button automatically detects a fall even if you can't press your button.*

- Comfortable & Lightweight
- Wearable as a pendant
- Waterproof, for use in places like the shower or pool

*Fall Button does not detect 100% of falls. If able, users should always push their help button when they need assistance. Fall Button is not intended to replace a caregiver. Equipment may vary from images shown. Prices subject to change without notice.

Only \$11/month!
(in addition to your plan)

Monitoring Center Phone Number

SAVE THIS NUMBER IN YOUR PHONE CONTACTS TO ENSURE YOU DON'T MISS AN IMPORTANT CALL FROM OUR MONITORING CENTER

1-800-932-3822

Please do not call
prior to reading
steps below.



Follow these steps below:

1. Add this phone number as a contact in your cell phone.
2. Ask your emergency contacts to do the same.
3. If you or your emergency contacts receive a call from this number, please answer!
4. If needed, call back the same number appearing on the Caller ID as not to cause any delay in processing an alarm.

Please note, there is no need to return the call unless requested by the monitoring center or if you have a question about an alarm.

SRDR-00026 Rev 2

Quick Start Guide

Easy steps to setup
your system today!



Complete Protection System (DUO)

We encourage our valued customers to use the testing processes outlined to allow operators more time to handle actual emergencies. Rest assured, a live operator will always be available to answer your call if an emergency were to occur.

SRDR-00026 Rev 2

QUICK SETUP

Your System is "ACTIVE" and READY TO USE


Complete Protection System

Includes the cellular base station, a Mobile Device, a neck pendant or wrist button and a cradle charger.




Plug the power cords into electrical outlets.


The electrical outlets should NOT be controlled by a light switch.



Place the Mobile Device on the Cradle Charger.



Turn on the ON/OFF switch on the back of the Base Station.



The **Emergency** and **Reset** buttons will light up in 10 seconds and the display screen will turn on. If the base station is set up correctly and connects to the cellular network, it will say "system ready" within 60 seconds.

You can use the volume control buttons (+ or -) on the front of the unit under the display screen if the voice is too loud or too soft.

Within a few minutes of turning on the system, the Cellular Base Station will obtain the date and time from the cellular network which will appear in the display area. You are not able to set the date and time manually.

If the Mobile Device is setup correctly and connects to the cellular network, it will also announce "system ready".


Next, be sure to test your system... see Testing page.

SRDR-00026 Rev 2

TESTING

Complete Protection System


1



Press and hold the test button

When the test button turns green, release it. This will place both devices in test mode.

2



Press the emergency button or pendant

The Base Station will announce, "Test Call sent to Emergency Response Center"
The Mobile Device will announce, "Call in Progress."

Always have your Mobile Device nearby when testing your system. It is possible that the Mobile Device will go to a live operator. If the operator speaks to you, please advise them that you are testing your system. If you do not explain to the operator that you are okay and that you are testing, EMS may be contacted.

3



If the test was successful...

The Base Station and Mobile Device will announce, "Thank you for testing your device"

4



If the test was unsuccessful...

The Base Station will announce, "User auto-test failed, please contact Customer Support"

PLEASE REFER TO THE TECHNICAL SUPPORT NUMBER ON THE INSERT THAT CAME WITH THIS GUIDE FOR ANY QUESTIONS OR CONCERNS.

If you have not set up your system... see Setup page first.

SRDR-00026 Rev 2